BLUE-2500-445

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DANGER LASER RADIATION!

Danger for eyes and skin

All protective measures for a safe operation of this laser must be applied.

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Safety instructions

Make sure to use the correct voltage

The operator has to make sure that laser radiation – also reflected laser radiation – higher than the maximum permissible level is avoided by technical or organisational measures. Always ensure, that the maximum permissible exposure (MPE) is not exceeded in the area accessible to the public.

In some countries, the operator must notify the accident insurance and the authority for industrial safety, before operating a laser. For more information, contact the relevant authorities.

Don't use the module if housing, or electrical parts are visibly damaged.

Do not turn device on and immediately off again frequently.

Always disconnect from mains, when handling, or installing the laser

Do not look into the laser beam directly, especially not with optical instruments.

Do not touch the device with wet hands.

When the laser diode becomes dim or broken, please contact your dealer timely.

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Laserworld Swisslas BLUE-2500-445 Lasermodul

High end Swisslas blue diode laser - used in collimation, scientific experimentation, optics instrument, laser display, and lighting etc.

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Data subject to change without notice. Errors and omissions expected!

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Connectors:

Caution: Do not switch supply voltage and ground – this will damage the driver!

You can supply TEC and Diode from one, or two separate PSUs – always connect both pins 3 and 5 to the PSU ground connection.


When using single ended modulation (one side connected to ground) connect pin2 to gnd!

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6 TEC 12V DC
5 TEC ground
4 Diode 12V DC
3 Diode ground
2 Modulation -
1 Modulation +

Mounting the module




For Swisslas laser modules the heat sink is mounted beneath the laser in the lower compartment of the housing. The images below show the mounting positions of a red/blue module in a Purelight unit.

Left image: The heat sink has to be directly beneath the laser module.

Use thermal grease under the TEC element of the laser (don't overturn the mounting screws, the TEC could break):

Centre image: Use 3 screws to mount the module (2 on one side, one centred on the other, see image), to avoid bending of the mounting plate. Use thin washers (0.1mm) to compensate for unevenness of the surface.

Right image: Use a good heat sink with a strong fan. Also use thermal grease between heat sink and mounting plate. The temperature of the module may not exceed 40°C – if it is exceeded the module will reduce output power (to ca. 50%) to protect the diodes!

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Laserworld (Switzerland) AG
Cleantech Europe GmbH